

POLYGLYCAN®-HV

**(Hyaluronic Acid + Chondroitin sulfates C4 & C6
+ N-acetyl-D-glucosamine Sterile Solution)**

**STERILE POST-SURGICAL LAVAGE, SYNOVIAL FLUID REPLACEMENT,
VISCOSUPPLEMENT**

**Veterinary Use Only
FOR USE IN HORSES**

Patent Nos.: 6,979,679 / 7,485,629 / 8,455,458



INDICATIONS

Polyglycan®-HV is indicated for the replacement of lost or damaged synovial fluid. Polyglycan®-HV contains naturally occurring components of the synovia that play a central part in maintaining the homeostatic environment of the joint.

BENEFITS

- **Helps replenish** lost synovial fluid with naturally occurring elements
- **Assists** in the normalization of synovial fluid viscosity
- **Post surgical lavage** of joint cavity and surfaces
- **Patented formulation**
- **Convenient pre-filled syringe**

CAUTION

This device is restricted to use by or on the order of a licensed veterinarian.

PACKAGING

LIST NO.	UNIT PACKAGE	CASE SIZE
1POL008	2.5 mL prefilled / 5mL syringe	40

Bimeda, Inc.

One Tower Lane, Suite 2250

Oakbrook Terrace, IL 60181 USA

Toll Free Tel. 888-524-6332 • Toll Free Fax. 877-888-7035

Email: Sales@BimedaUS.com

See reverse side for additional information



BimedaUS.com

POLYGLYCAN®-HV

*(Hyaluronic Acid + Chondroitin sulfates C4 & C6 +
N-acetyl-D-glucosamine Sterile Solution)*

**STERILE POST-SURGICAL LAVAGE, SYNOVIAL FLUID REPLACEMENT,
VISCOSUPPLEMENT**

Veterinary Use Only
FOR USE IN HORSES

Patent Nos.: 6,979,679 / 7,485,629 / 8,455,458

CAUTION:

This device is restricted to use by or on the order of a licensed veterinarian.

DESCRIPTION:

Polyglycan®-HV is a patented formulation designed to replace lost or damaged synovial fluid. Polyglycan®-HV contains naturally occurring components of the synovia that play a central role in maintaining the homeostatic environment of the joint.

COMPOSITION:

Highly viscous aqueous solution of defined fractions of purified Hyaluronic acid, Chondroitin sulfates C4 & C6 in a 5% solution of N-acetyl-D-glucosamine.

EACH VIAL CONTAINS:

Hyaluronic acid sodium salt.....	20 mg
Sodium chondroitin sulfate.....	250 mg
N-acetyl-D-glucosamine.....	125 mg

STORAGE:

Refrigerate whenever possible. DO NOT FREEZE. Maintain product at or below 65°F (18°C) during shipment and storage.

PROPERTIES:

Glycosaminoglycans are important components of all extracellular tissue structures including cartilage and synovial fluid. The active components in Polyglycan®-HV exhibit viscoelastic and polyionic properties similar in nature to synovial fluid.

Introduce Polyglycan®-HV into the synovial space as a post-surgical lavage, synovial fluid replacement or viscosupplement to assist in the normalization of the joint.

DIRECTIONS:

Instill Polyglycan®-HV into the synovial space to replenish lost synovial fluid. Quantity to be determined by Veterinarian Surgeon. Discard any unused portion of vial contents.

WARNING:

Do not administer to animals that are to be slaughtered for use in food. Keep out of reach of children.

